

Wystopper- a flexible clear soft stoma stopper designed by a clinician to maintain stoma sites, especially continent urinary diversions and antegrade continence enema (ACE) tracts.




### Tapered Tip

The clear soft silicone aids patient comfort and tapered tip to ease insertion.



### Mushroom Cap

The mushroom shaped cap lies flush against the skin for discreetness.



### Choice of sizes

Available in sizes 10Ch, 12Ch, 14Ch, and 16Ch in 30mm or 60mm lengths.

DESCRIPTION	10Ch	12Ch	14Ch	16Ch	QTY
Wystopper Stoma Stopper 30mm	WS0310	WS0312	WS0314	WS0316	3
Wystopper Stoma Stopper 60mm	WS0610	WS0612	WS0614	WS0616	3